

# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

For additional information, please visit <https://www.woundcarehandbook.com/>

For Greenwich Integrated Woundcare Pathway, please refer to: [Integrated wound care pathway](#)

Prices listed as per June 2025

## 1. PRIMARY DRESSINGS (Require a secondary dressing)

**Lightly exuding or granulating wounds. Low adherence to wound bed. Non-absorbent**

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
Atrauman®	7.5 x 10	£0.42	10	5	369-7547	Non-adherent polyester mesh wound contact layer. Suitable for fragile skin that is likely to tear easily. Effective for up to 7 days Petrolatum free Impregnated with triglycerides	
	10 x 20	£0.97	10	5	281-3046		
Mepitel One®	9 x 10	£2.52	5	5	350-6896	If clinically indicated up to 14 days	
<b>Alginate dressing</b>							
Kaltostat®	5 x 5	£1.05	10	5	007-0789	Kaltostat® is haemostatic and useful for bleeding wounds. Remove when bleeding stops. If clinically indicated up to 7 days	

## 2. ABSORBENT PERFORATED DRESSING WITH ADHESIVE BORDER

**Low to moderate exuding wounds, protective covering for primarily closed and postoperative wounds. Low adherence primary or secondary dressing.**

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
Cosmopor E®	5 x 7.2	£0.09	50	10	374-0115	Self-adhesive island wound dressing with non-adherent absorbent pad. If clinically indicated – up to 7 days	NOT showerproof
	8 x 10	£0.19	25	10	374-0123		
	8 x 15	£0.31	25	10	3740131		
	10 x 20	£0.50	25	10	374-0149		

# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

For additional information, please visit <https://www.woundcarehandbook.com/>

For Greenwich Integrated Woundcare Pathway, please refer to: [Integrated wound care pathway](#)

Prices listed as per June 2025

## 3. HYDROCOLLOID DRESSING

All tissue types, but only those with low to moderate exudate. Will rehydrate necrosis and slough. Promote granulation. Only use as a primary dressing.

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
DUODERM EXTRA THIN®	10 x 10	£1.45	10	5	027-7897	Only prescribe a small amount as the wound exudate will increase as the wound debrides and hydrocolloids will no longer be suitable. <b>Notes for Care home –</b> <b>DO NOT apply to diabetic/ischaemic foot ulcers</b> <b>Avoid in infected wounds</b> <b>Up to 7 days use</b>	<ul style="list-style-type: none"> <li>Do not allow to dry out.</li> <li>Do not pack into sinuses as risk of plug formation.</li> <li>Dry necrotic wounds 3rd degree burns.</li> <li>Heavily bleeding wounds. Caution in Fungating tumours</li> </ul>

## 4. SIMPLE FOAM DRESSINGS

Granulating, epithelising or sloughy wounds with moderate to heavy exudate. Can be primary or secondary.

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
AQUACEL FOAM PRO ADHESIVE®	8 x 8 10 x 10 15 x 15	£1.54 £2.36 £4.24	10 10 10	5 5 5	422 359 422 357 422 358	For exuding and non-exuding wounds. Do not cover foams with occlusive dressings. May be used for wounds such as leg, pressure and diabetic foot ulcers, surgical wounds, traumatic wounds and partial thickness burns If clinically indicated – up to 7 days use	Not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or components Not compatible with oxidising agents such as hydrogen peroxide or hypochlorite solutions.
AQUACEL FOAM non-adhesive	10 x 10 15 x 15	£2.78 £4.67	10 10	5 5	369-9329 370-2750	For moderate to highly exuding wounds Can be used with secondary dressings or compression bandages If clinically indicated – up to 7 days use Care Homes might consider using Activheal foam if more suitable.	

# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

For additional information, please visit <https://www.woundcarehandbook.com/>

For Greenwich Integrated Woundcare Pathway, please refer to: [Integrated wound care pathway](#)

Prices listed as per June 2025

<b>KLINIDERM FOAM® Silicone non- adhesive</b>	10 x 20	£1.74	10	5	394-7223	Absorbent foam dressing with a soft silicone wound contact layer which protects the surrounding skin and reduces pain during dressing changes. The absorbent foam layer absorbs and retains exudate. In addition, wound exudate partly evaporates through the dressing, increasing the total absorption capacity. Supports moist wound treatment. If clinically indicated up to 7 days	
<b>Tegaderm® Foam adhesive</b>	10 x 11 14.3 x 15.6	£2.48 £4.57	10 5	5 5	304-1589 304-1605	Tegaderm may be preferred for wounds in 'difficult to dress' areas. If clinically indicated up to 7 days	Not to be used on third-degree burns, dry or non-draining wounds, and as a pressure reduction device or surgical sponge. Avoid using with oxidising agents such as hydrogen peroxide.
<b>Urgostart® Plus pad</b>	10 x 10	£5.14	10	5	406-4440	Dressing may be changed every 2 to 4 days and left in place for up to 7 days depending on the level of exudate and the condition of the wound. The recommended treatment duration is a minimum of 4 to 5 weeks.	Not to be used in mucous membrane wounds, malignant wounds, fistula wounds, epidermolysis bullosa wounds, necrotic wounds and haemorrhagic wounds.
<b>Urgostart® Plus border</b>	13 x 13	£8.84	10	5	406-4416	To be used when patients are in compression hosiery Generally changed every 1 to 2 days during the wound desloughing phase and then adapted based on the wound's exudate level and clinical progress, but can be left in place for up to 7 days.	Not to be used in heavily bleeding wounds, cancerous (malignant) wounds, and wounds with deep abscesses. Not recommended for epidermolysis bullosa, dry necrosis, or known sensitivity to the dressing components

# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

For additional information, please visit <https://www.woundcarehandbook.com/>

For Greenwich Integrated Woundcare Pathway, please refer to: [Integrated wound care pathway](#)

Prices listed as per June 2025

## 5. SUPER ABSORBENT DRESSINGS

Management of heavy exuding wounds. Can be primary or secondary dressing.

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
Kerramax® Care Border Adhesive	16 x 16	£4.70	10	5	379-1977	For heavily exuding wounds, can increase wear time. If clinically indicated up to 7 days use.	Must not be cut or torn The high absorption capacity can lead to dehydration when used in areas where moisture is essential
	16 x 26	£7.62	10	5	379-2124		

## 6. ANTIMICROBIAL DRESSINGS

### 6.1 Enzyme Alginogel

For infected wounds – Short-term use 7 - 14 days usually advised until infection resolves. Require a secondary dressing.

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
For moderate to highly exuding wounds							
Flaminal® Forte	15g	£8.58	1	1	324-2963	Ensure wound bed is covered with 4-5mm of Flaminal® gel	Not be used on full-thickness burns (third-degree burns) or known sensitivity to alginate dressings or polyethylene glycol
	50g	£28.42	1	1	344-9592		
Flaminal® Hydro	15g	£8.44	1	1	324-2971		
	50g	£28.42	1	1	344-9600		

### 6.2 IODINE DRESSINGS

For infected wounds – maximum single application 50g or maximum weekly application 150g. All require a secondary dressing.

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
For Superficial/Low exuding wounds							
INADINE®	9.5 x 9.5	£0.51	10	10	037-1229	Effective against MRSA but not pseudomonas	• DO NOT USE if: -Sensitivity to Iodine

# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

For additional information, please visit <https://www.woundcarehandbook.com/>

For Greenwich Integrated Woundcare Pathway, please refer to: [Integrated wound care pathway](#)

Prices listed as per June 2025

						If wound not improving to treatment after 4 weeks treatment refer to Tissue Viability Service If clinically indicated – up to 7 days	<ul style="list-style-type: none"> <li>-Before or after use of radio-iodine</li> <li>-Children, pregnant or lactating women</li> <li>-Renal</li> <li>• Use under medical supervisions: <ul style="list-style-type: none"> <li>-Thyroid disorders</li> <li>-Diabetes taking sulphonylureas</li> <li>-Lithium</li> </ul> </li> </ul>
<b>For Heavily exuding/necrotic wounds</b>							
Iodoflex <sup>®</sup> paste	5g	£4.54	5	5	073-1547	<p><b>ONLY to be used under Tissue Viability Service supervision</b></p> <p>Adherent crust may form when the dressing is changed with significant frequency.</p>	<ul style="list-style-type: none"> <li>• DO NOT USE if: <ul style="list-style-type: none"> <li>-Sensitivity to Iodine</li> <li>-Hashimoto's thyroiditis</li> <li>-non-toxic nodular goitre children</li> <li>- Pregnant or lactating women</li> </ul> </li> <li>• Use under medical supervisions: <ul style="list-style-type: none"> <li>-Thyroid disorders</li> <li>-Diabetes taking sulphonylureas</li> <li>-Lithium</li> </ul> </li> </ul>

## 6.3 SILVER DRESSINGS

**For infected wounds – Short-term use 7 – 14 days usually advised until infection resolves. All require a secondary dressing.**

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Max of days dressings can stay on	Cautions/ Contraindications
Aquacel Ag + Extra <sup>®</sup>	5 x 5 10 x 10	£2.15 £5.11	10 10	5 5	386-2703 386-2695	<ul style="list-style-type: none"> <li>• 2 weeks use only</li> <li>• Effective up to 7 days</li> <li>• Can be used in combination with traditional swabs, foam</li> </ul>	Dressings should be changed every 3-7 days If the wound needs daily change, do not	<ul style="list-style-type: none"> <li>• DO NOT USE if known sensitivity to silver</li> <li>• Must be removed prior to x-ray, MRI, ultrasound or diathermy</li> </ul>
Urgotul <sup>®</sup> Silver	10 x 12	£3.39	16	5	329-5995			

# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

For additional information, please visit <https://www.woundcarehandbook.com/>

For Greenwich Integrated Woundcare Pathway, please refer to: [Integrated wound care pathway](#)

Prices listed as per June 2025

<b>Aquacel Ag + Ribbon®</b>	1 x 45	£3.37	5	5	386-2729	dressings and absorbent dressings • Urgotul Silver is cost-effective • Aquacel Ag Ribbon – primarily for cavity/sinus wounds with moderate to heavy exudate	use Silver dressings. Use standard dressings and consider systemic antibiotic if indicated or leave inner silver dressing intact and change soiled outer dressing.	
<b>Acticoat® flex 3</b>	10 x 10	£9.61	5	5	351-0799			

## 7. ODOUR ABSORBING DRESSINGS

Non-adhesive secondary dressing designed to eliminate unpleasant odour

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
<b>Cinesteam®</b>	11 x 19	£6.33	1	1	423-0447	Consisting of an absorbent layer, together with a sachet containing cinnamon - offers a triple action anti-odour approach Absorbs excess exudate If clinically indicated, up to 7 days use.	Not be used on patients with known sensitivities to any of its components, especially cinnamon

## 8. FILM DRESSINGS

Primary/secondary vapour permeable film for shallow, dry, non-infected wounds e.g. donor sites, minor burns, grazes, postoperative incisions.

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	Additional information	Cautions/ Contraindications
<b>365 film dressings®</b>	6 x 7 10 x 12	£0.07 £0.16	100 50	10 10	Showerproof.	-

# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

For additional information, please visit <https://www.woundcarehandbook.com/>

For Greenwich Integrated Woundcare Pathway, please refer to: [Integrated wound care pathway](#)

Prices listed as per June 2025

					Stretch opposite corners for easy removal. Dressings can be left in place for up to 7 days, if signs of infection are observed, remove the dressing.	
--	--	--	--	--	---	--

## 9. PHYSICAL DEBRIDEMENT PADS

Debridement of superficial wounds containing loose slough and debris, and for removal of hyperkeratosis from the wound.

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
Alprep® Pad	7 x 9	£6.78	5	5	417-1146	Used for manual debridement of slough/hyperkeratosis	Not to be used on wounds with suspected or active osteomyelitis, infected wounds requiring topical or systemic antibiotics, or wounds with 100% healthy granulation tissue

## 10. SKIN PREPARATION

Skin barrier and moisturiser for intact and injured skin.

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
Dressit Dressing Pack	S/M M/L	£0.70	1	1 pack	324-3961 301-0675	Only to be used for sterile wound care When considering prescribing/recommending a dressing pack, remember:  Dressing packs should only be used where clinically indicated and in line with employing organisation policies, procedures and guidelines	There is no component of any dressing pack which should be left in permanent contact with the wound Dressing packs containing cotton wool balls should not be used due to fibre-shed.

# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

For additional information, please visit <https://www.woundcarehandbook.com/>

For Greenwich Integrated Woundcare Pathway, please refer to: [Integrated wound care pathway](#)

Prices listed as per June 2025

<b>Clinimed® Guaze</b>	10 x 10	£0.88	100	1 pack	-	For cleansing of chronic wounds, where a sterile wound pack is not required	
<b>Derma Protective® Plus Skin Protectant</b>	115g	£4.99	1	1	399-6931	Pea to Palm sized amount Do NOT over apply & ensure the skin is visible through the cream. For moisture associated skin damage <b>NOT</b> to relieve pressure.	
<b>Esenta® sting free skin barrier pump spray</b>	28ml	£6.17	1	1	423286	Sting-Free Skin Barrier creates a breathable barrier on skin that protects skin for up to 72 hours. Reduces the risk of skin damage and irritation caused by body waste, enzymes, and adhesives, and is safe and effective with frequent use.	
<b>Lifteez® non- sting medical adhesive remover</b>	50ml	£5.60	1	1	389-7147	Removal of dressing or device from fragile skin. Spray on porous dressing for removal. Non-Porous - lift edge of dressing with gradual spray underneath	
<b>Epimax® original cream</b>	500g	£2.72	1	1	393-4270	First choice as a soap substitute or emollient	The Medicines and Healthcare products Regulatory Agency (MHRA) has identified a significant fire hazard risk associated with emollients and dressings, particularly when they are applied to clothing or bedding and then dried. See <a href="#">MHRA advice</a> and <a href="#">SEL emollient guide</a> for more information
<b>Do not use Conotrane® as barrier cream.</b>							



# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

For additional information, please visit <https://www.woundcarehandbook.com/>

For Greenwich Integrated Woundcare Pathway, please refer to: [Integrated wound care pathway](#)

Prices listed as per June 2025

## 11. ABSORBENT AND SUPERABSORBENT PADS

For heavily exuding wounds, such as leg ulcers, pressure ulcers, and diabetic foot ulcers. Also work well with compression therapy

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
<b>DRYMAX® Super</b>	10 x 11	£0.91	10	5	421-3369	More absorbent than Kliniderm® superabsorbent. For more highly exuding wounds.	Absorb excess fluid from wounds, so should not be used on wounds that are already dry. As the dressing swells when it absorbs fluid, it should not be used in areas where swelling could be problematic, such as wound cavities.
	11 x 20	£1.09	10	5	421-3377		
	20 x 20	£1.93	10	5	421-3385		
	20 x 30	£2.44	10	5	421-3393		
	37 x 56	£5.02	10	5	421-3401		

## 12. BANDAGES

### 12.1 TUBULAR RETAINING BANDAGE

Dressing retention if retention bandage is unsuitable

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
<b>Tubifast® Blue</b>	7.5 x 3m	£2.26	1	5	327-6433		
	7.5 x 5m	£3.77	1	5	327-6441		
<b>Tubifast® Yellow</b>	10.75 x 3m	£3.71	1	5	327-6466		
	10.75 x 5m	£6.09	1	5	327-6474		

# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

For additional information, please visit <https://www.woundcarehandbook.com/>

For Greenwich Integrated Woundcare Pathway, please refer to: [Integrated wound care pathway](#)

Prices listed as per June 2025

## 12.2 SUB-COMPRESSION WADDING BANDAGE

Soft absorbent protection of skin under compression bandaging and for padding and reshaping the lower limb as required.

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
K-SOFT®	10 x 3.5m	£0.49	1	5		Part of K-FOUR multilayer compression bandage system	Not recommended for patients with severe arterial disease (ABPI <0.6), diabetic microangiopathy, ischaemic phlebitis, septic thrombosis, and ulceration caused by infection.

## 12.3 LIGHT SUPPORTING BANDAGE

Light compression bandage for venous leg ulcers

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
K-Plus®	10 x 8.7	£2.47	1	5	248-1232	Type 3a light compression bandage that is a highly extensible elastic bandage and donates up to 20mmHg of sub-bandage pressure at the ankle. Third layer in the K-Four multilayer compression system for venous leg ulcers.	Not recommended for patients with severe arterial disease (ABPI <0.6), diabetic microangiopathy, ischaemic phlebitis, and septic thrombosis

## 12.4 LIGHT SUPPORTING CONFORMING BANDAGE

Reducing mild oedema, non-venous leg ulcers or where compression cannot be tolerated.

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
K-LITE®	10 x 4.5	£1.09	1	5		Part of K-FOUR multilayer compression bandage system	Known hypersensitivity to nylon, viscose, elastomer. Not recommended for patients with arterial disease (ABPI < 0.8), diabetic microangiopathy, ischaemic phlebitis, septic thrombosis

# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

For additional information, please visit <https://www.woundcarehandbook.com/>

For Greenwich Integrated Woundcare Pathway, please refer to: [Integrated wound care pathway](#)

Prices listed as per June 2025

K-SOFT® Long	10 x 4.5	£0.62	1	5			Not recommended for patients with severe peripheral arterial disease (ABPI < 0.6), diabetic microangiopathy, or acute deep vein thrombosis
--------------	----------	-------	---	---	--	--	--

## 12.5 ZINC PASTE BANDAGE

For the treatment of venous ulcers.

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
ZipZoc Zinc® oxide stockings	1 stocking	£15.90	4	1	343-3739	Should be changed every 3 to 5 days Patch test before use	Hypersensitivity to any ingredient of the stocking and acute eczematous lesions
Viscopaste® PB7 (10%)	7.5 x 6m	£3.99	1	1	033-2734	Patch test before use	Hypersensitivity to any ingredient in the bandage.

## 13. COMPRESSION SYSTEMS

### 13.1 TWO LAYER SYSTEMS – (NOT TO BE APPLIED IN CARE HOMES)

Compression therapy for venous leg ulcers with ABPI of > 0.8 and venous oedema

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
UrgoKTwo® KIT (40mmHg)	18-25 x 10m 25-32 x 10m	£8.72 £9.53	1 1	1 1	043-5891 327-4685	• MUST only be applied by competent practitioners • Suitable for ankle circumference of < 25 cm	• DO NOT USE: - Hypersensitivity to latex - ABPI < 0.8 - Arterial disease
UrgoKTwo® REDUCED KIT (22mmHg)	18-25 x 10m 25-32 x 10m	£8.72 £9.53	1 1	1 1	333-8480 360-2869	• Can be washed and reused • Latex-free version available • Can be left for up to 7 days	• Use under medical supervision for: Diabetes, heart failure or advanced microvascular disease patients

# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

For additional information, please visit <https://www.woundcarehandbook.com/>

For Greenwich Integrated Woundcare Pathway, please refer to: [Integrated wound care pathway](#)

Prices listed as per June 2025

## 13.2 TWO LAYER SHORT STRETCH – (NOT TO BE APPLIED IN CARE HOMES)

Compression therapy of limbs with venous disorders (i.e. venous leg ulcers) and the control and reduction of chronic oedema and lymphoedema in the lower limbs and arms. This form of compression therapy can be used with caution in patients presenting with mixed aetiology leg ulcers, as the wound may not be a pure venous ulcer but may be mixed with, for example, arterial disease.

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
ACTICO®	8 x 6m	£3.54	1	1	314-0886	<ul style="list-style-type: none"> <li>• MUST only be applied by competent practitioners</li> <li>• Confirm ABPI doppler reading before applying</li> <li>• Contains low sensitivity latex</li> <li>• Can be left for up to 7 days</li> </ul>	<ul style="list-style-type: none"> <li>• DO NOT USE: <ul style="list-style-type: none"> <li>- Hypersensitivity to latex</li> <li>- ABPI &lt; 0.8 - Arterial disease</li> </ul> </li> <li>• Use under medical supervision for: Diabetes, heart failure or advanced microvascular disease patients</li> </ul>
	10 x 6m	£3.67	1	1	271-5431		
	12 x 6m	£4.68	1	1	314-0894		

## 13.3 COMPRESSION WRAP

For keep dressings, devices and creams in place without causing discomfort to the patient

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	Additional information	Cautions/ Contraindications
L&R ReadyWrap®	S,M,L,XL,XXL	£16.01	1	1	Adjustable, short-stretch compression garments with VELCRO™ brand fasteners are designed for an easy-to-use compression therapy for patients with Venous Leg Ulcers, Venous disease, Chronic oedema and Lymphoedema. Available in Black or Beige 40-50mmHg compression	<p>The system may be unsuitable when there is heavy exudate and for large wounds. Caution is required when:</p> <ul style="list-style-type: none"> <li>• Cardiac overload is suspected</li> <li>• Patients have diabetes</li> <li>• Patients have advanced small vessel disease</li> <li>• Arterial disease is present</li> <li>• Renal failure is present</li> <li>• Rheumatoid arthritis is present</li> <li>• Following application if the patient</li> </ul>

# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

For additional information, please visit <https://www.woundcarehandbook.com/>

For Greenwich Integrated Woundcare Pathway, please refer to: [Integrated wound care pathway](#)

Prices listed as per June 2025

						shows signs of pain reduced colour or numbness to their extremities then it should be removed immediately Do not use if the ulcer is below the malleolus bone.
--	--	--	--	--	--	---

## 14. COMPRESSION THERAPY - if any queries from lymphoedema team, contact [oxl-tr.greenwichlymph@nhs.net](mailto:oxl-tr.greenwichlymph@nhs.net)

### Treatment for venous disorders, chronic oedema and lymphoedema

**Education on how to wear and wash the stockings are essential before prescribing.**

Measurements of widest part of the calf, ankle and foot length required before prescribing. (Videos for patient: [medi UK how to measure for mediven® off the shelf](#))

Both clinicians and patients can access the information on measurement:

[Measurement form for Mediven](#), [Measurement form for Duomed](#), [Help With Measuring \(jobst.co.uk\)](#)

Ensure patients/carer can use these products before prescribing. **Issue 1 pack initially for trial on each affected limb**

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	Additional information	Cautions/ Contraindications
<b>Mediven Plus®</b> <b>Class 1*</b> (18-21mmHg) <b>Class 2*</b> (23-32mmHg)	(Open/ closed toe) (Black/beige) Size 1 -7	£31.52 £31.52	2 2	2 pair 2 pair	Can last for up to 6 months For telangiectasia or reticular veins, oedema, prevention of ulceration or recurrence (Petite and extra wide options available) Used for chronic venous disease stages <a href="#">Cos-c3</a> RAL compression classification	Class 1 for all conditions Class 2 not recommended for patients with arterial disease (ABPI< 0.8)
<b>Mediven for men®</b> <b>Class 1</b> (18-21mmHg) <b>Class 2</b> (23-32mmHg)	(Closed toe only) (Navy/Grey/black) Size 1 -7	£31.06 £31.41	2 2	1 pack 1 pack	Can last for up to 6 months For Broader feet and early to moderate venous disease, soft oedema (Petite options available) Used for chronic venous disease stages <a href="#">Cos-c3</a> RAL compression classification	

# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

For additional information, please visit <https://www.woundcarehandbook.com/>

For Greenwich Integrated Woundcare Pathway, please refer to: [Integrated wound care pathway](#)

Prices listed as per June 2025

<b>Mediven Mondi® MTM Class 1 (18-21mmHg) Class 2 (23-32mmHg)</b>	<p>Made to measure</p>	<p>£37.75</p> <p>£39.53</p>	<p>2</p> <p>2</p>	<p>2 pairs (4 garments)</p> <p>2 pairs (4 garments)</p>	<p>Prescription for 2 garments for each leg. Measurements to be added to a <a href="#">custom-fit form</a>. Flat-knitted medical compression stocking used for compression of the lower extremities, mainly for the treatment of disorders of the lymphatic system:</p> <ul style="list-style-type: none"> <li>• Primary and secondary lymphoedema stages 1 and 2</li> <li>• Lipoedema stages 1 und 2</li> <li>• Post-operative oedema</li> <li>• Post-traumatic oedema</li> <li>• Scar treatment</li> </ul> <p>6 months usage RAL compression classification</p>	<p>Advanced peripheral arterial occlusive disease, decompensated heart failure septic thrombophlebitis, and phlegmasia cerulea dolens. Additionally, any local leg condition where compression would cause further damage (fragile skin, dermatitis, or gangrene).</p>
<b>Duomed soft® Class 1 (14-17mmHg) Class 2 (18-24mmHg)</b>	<p>(Sand/Black) S-XXL</p>	<p>£22.82</p> <p>£22.82</p>	<p>2</p> <p>2</p>	<p>1 pair</p> <p>1 pair</p>	<p>Can last for up to 3 months Contains a softer band at the top preventing indentation and soreness below the knee British Standard compression classification</p>	<p>British Standard Compression: Class 1 (14-17mmHg) for all conditions Class 2 (18-24 mmHg) not recommended for arterial disease (ABPI&lt; 0.8)</p>
<b>Jobst UlcerCare®</b>	<p>From Small to XXXX Large</p>	<p>£34.29</p>	<p>2 liners 1 outer</p>	<p>1 pack</p>	<p>Can last up to 3-6 months. JOBST UlcerCare® is a 2-in-1 compression system designed to effectively manage venous leg ulcers. The two-part system provides graduated compression with a combined nominal compression of 40 mmHg. RAL compression classification</p>	
<b>Jobst UlcerCare® Compression liner pack</b>	<p>From Small to XXXX Large</p>	<p>£20.71</p>	<p>3</p>	<p>1 pack</p>	<p>Can last for up to 3 months Liners only</p>	<p>Does not contain compression stockings.</p>

# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

For additional information, please visit <https://www.woundcarehandbook.com/>

For Greenwich Integrated Woundcare Pathway, please refer to: [Integrated wound care pathway](#)

Prices listed as per June 2025

<b>Jobst Elvarex®</b> <b>Made to measure of Custom fit</b> <b>Class 1</b> (18-21mmHg) <b>Class 2</b> (23-32mmHg)	Custom fit knee high	£37.52	1	2	Can last up to 6 months Have a higher static stiffness and a higher interservice pressure. Custom-fit to the anatomical shape of the limb when patients do not fit into the ready to wear sizing RAL compression classification	Severe arterial insufficiency, uncontrolled congestive heart failure, and untreated infection
		£39.29	1	2		
<b>Jobst Elvarex®</b> <b>Soft Made to measure of custom fit</b> <b>Class 1</b> (18-21mmHg) <b>Class 2</b> (23-32mmHg)	Custom fit knee high	£36.23	1	2		
		£38.50	1	2		
<b>Ulcer X® Kit</b>	Small to X Large	£26.98	-	1 pack	Can last for up to 6 months Contains of 1 stocking and 2 liners Two part system providing 45mmHg of graduated compression US compression level	Contraindicated when: • Ischemia (e.g. advanced arterial disease) • Untreated septic phlebitis • Uncontrolled congestive heart failure • Phlegmasia cerulea dolens
<b>Ulcer X Refill kit</b>	Small to X Large Closed toe only	£26.55	4	1 pack	Protection: due to the patented liner, the wound dressing also remains fixed in place at night, and the wound is not affected adversely when used correctly -neither when donning or doffing Effective: better healing rates in small and medium-sized ulcers compared to bandages The liner offers 15-20mmHg of compression	



# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

<b>Sigvaris® Traditional 500 Class 2</b> (23-32mmHg)	(Open toe only) Small to X Large	£30.21	2	1 pack	Can be used up to 6 months. Suitable for those who need consistent compression The unique stress-strain behaviour of the yarn composition ensures high pressure stability throughout the entire day and effective treatment even in advanced venous diseases. Easy to do even at high compression levels RAL compression classification	Contraindicated in: <ul style="list-style-type: none"> <li>Severe peripheral arterial disease, including low ABPI, those who have undergone a peripheral arterial bypass graft, and those with significant peripheral neuropathy.</li> <li>uncontrolled congestive heart failure (CHF)</li> <li>acute dermatitis</li> <li>signs of infection in the legs.</li> </ul>
<b>Sigvaris® Essential comfortable Class 1</b> (18-21mmHg)	(Closed toe/open toe) Small to X Large	£30.08	2	1 pack	Can be used up to 6 months Suitable for those with sensitive skin and for controlling skin temperature. RAL compression classification	Contraindicated in: <ul style="list-style-type: none"> <li>severe peripheral arterial occlusive disease</li> <li>advanced neuropathy, or significant sensory disturbances in the extremity.</li> <li>Avoid in certain inflammatory conditions or skin diseases with massive exudation</li> </ul>
<b>Haddenham Doktus® MTO Class 2</b> (23-32mmHg) <b>Class 3</b> (34-46mmHg)	Made to order	£30.17 (both classes)	1	2	Can be used for 6 months Available in four colours, with grip top options, closed toe, open toe and footless. Short, regular, and long foot and length can be extra wide RAL compression classification	Class 2 - Severe peripheral arterial disease, acute deep vein thrombosis, and certain infections or wounds. Class 3 - Severe arterial insufficiency, untreated or active DVT, severe peripheral neuropathy, significant cardiac issues and certain skin conditions.
<b>Haddenham star cotton Class 1</b> (18-21mmHg) <b>Class 2</b> (23-32mmHg)	Short, regular, long	£29.85 (both classes)	1	1 pair (2 garments)	Can last up to 6 months RAL compression classification	Class 1 - Advanced arterial disease in the legs, severe heart disease, and untreated septic phlebitis Class 2 - severe arterial insufficiency, uncontrolled congestive heart failure, and untreated infection



# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

For additional information, please visit <https://www.woundcarehandbook.com/>

For Greenwich Integrated Woundcare Pathway, please refer to: [Integrated wound care pathway](#)

Prices listed as per Junel 2025

Further information:

Mediven - [Prescriber Support](#)

Mediven – how to take care of compression stockings - [Taking care of compression stockings – how do you wash them correctly? | medi](#)

Sigvaris - [How to measure for compression products | Sigvaris.com](#)

Haddenham - [Haddenham Custom Garments Measuring Guide.pdf](#)

## 15. TAPES

### Fixation of dressings, swabs, catheters and tubes

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
Clinipore®	2.5 x 5m	£0.52	5	5	299-0109		-

## 16. IRRIGATION AND CLEANSING

### Topical irrigation and cleansing of traumatic and surgical wounds

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	Additional information	Cautions/ Contraindications
(First line) WARM TAP WATER	-	£0.00	-	-	NICE GUIDELINE: <a href="#">Surgical site infections: prevention and treatment</a>	The decision to use tap water to cleanse wounds should take into account the quality of water, nature of wounds and the patient's general condition, including the presence of comorbid conditions

# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

For additional information, please visit <https://www.woundcarehandbook.com/>

For Greenwich Integrated Woundcare Pathway, please refer to: [Integrated wound care pathway](#)

Prices listed as per June 2025

## 17. PODIATRY USE ONLY ITEMS

### Treatment of diabetic foot ulcers and venous leg ulcers as per NICE TA [MTG 42]

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
URGOSTART CONTACT®	10 x 10	£4.45	10	5	386-1390	<b>For podiatry use ONLY</b> Should be changed every 1 to 2 days during the desloughing phase, then as often as required depending on the exudate volume and the clinical wound condition	-

## 18. OTHERS

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	PIP code	Additional information	Cautions/ Contraindications
AQUACEL Extra®	5 x 5cm	£1.09	10	5	ELY377	Primary absorbent, high tensile strength for cavity and sinus wounds For moderate to highly exudating wounds – requires a second dressing	
	10 x 10cm	£2.59	10	5	ELY378		
	15 x 15cm	£4.87	10	5	ELY379		
AQUACEL RIBBON®	2 x 45	£2.69	5	2	ELY013	Hydrofibre - primary dressing for cavity/sinus wounds with moderate heavy exudate	Should not be packed too tightly – for cavity wounds only
KerraPro® sheet	10 x 10 x 1.2	£13.67	1	1	3756061	Patients at risk of developing pressure damage on bony prominences or on areas where medical devices may cause damage to the skin due to pressure. Should be used on intact skin on areas where damage is likely to occur or where a category I pressure ulcer is already present. Skin underneath KerraPro® should be inspected regularly. Can be washed with soap and water, dried thoroughly and reused on the same patient.	Known sensitivity to silicone Should not be used as a wound dressing or directly on broken skin.

# Greenwich Combined Wound Care Formulary

(Additional information and Max number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing)

For additional information, please visit <https://www.woundcarehandbook.com/>

For Greenwich Integrated Woundcare Pathway, please refer to: [Integrated wound care pathway](#)

Prices listed as per June 2025

## 19. SCAR TREATMENT PREPARATIONS (PRESCRIBE ONLY UNDER SECONDARY CARE RECOMMENDATIONS)

For keloid and hypertrophic scars less than 2 years old resulting from trauma, burns or surgery

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Max Qty on Rx	Additional information	Max of days dressings can stay on	Cautions/ Contraindications
ADVASIL CONFORM®	10 x 10	£6.04	1	1	Prescribe one gel sheet only, silicone gel sheets can be washed and reused.	Once the gel sheet becomes difficult to clean, a new piece should be used. This usually occurs after 14 to 28 day intervals.	DO NOT USE: •Hypersensitivity to silicone •Broken skin •Open wounds •Any dermatological conditions e.g. eczema, acne
SCARSIL® GEL	30g tube	£16.92	30g	30g	For scars on the face, ears and neck are highly visible	-	

## 20. SELF CARE ITEMS (DO NOT PRESCRIBE)

Proprietary Name	Sizes (cm)	Unit cost	Pack size	Additional information	Cautions/ Contraindications
Seal Tight	All types and sizes	OTC	-	Physical cover to enable user to shower or bathe when using a cast or dressing for broken bones and wounds. Patients are to <b>buy over the counter</b>	No known contraindications
Limb-O® waterproof protectors	All types and sizes	OTC	-		
Kerraped Plus boots	All sizes	OTC	1	<u>DUAL STATUS:</u> For offloading pressure from the front of the foot – <b>to buy over the counter</b> Neuropathic foot ulcers or patients with bulky bandages who cannot wear their own footwear refer to <a href="#">Specialist-Formulary (selondonccg.nhs.uk)</a>	No known contraindications
Benefoot® Medical Shoe – Original post op	Men M/L Women S/M/L	OTC	1	To aid compliance with compression – <b>to buy over counter</b>	Not to be used for housebound patients

## Greenwich Prescription Request from Care Home Nurses for selected dressings from the Primary Care Wound Formulary

**To be prescribed by the GP**

GP/ Surgery:	Site: (left/right)
HCP requesting prescription: HCP name in print:	Frequency of change:        /week
Date requested:	Type of wound: light / medium / heavy exudate /other <b>please be specific</b>
Patient name:	Assessed by (HCP name in print):
DOB:	Next review date:

Dressing name	Size (delete unwanted size)	Quantity of dressing	Comment/alternative size	Additional information
DRYMAX® Super	10 x 11cm	5		If clinically indicated - up to 7 days use
	11 x 20cm	5		
	20 x 20cm	5		
	20 x 30cm	5		
	37 x 56cm	5		
Kliniderm® Foam Silicone non-adhesive	10 x 20cm	5		If clinically indicated – up to 7 days use
Aquacel Foam Pro Adhesive®	8 x 8cm	5		If clinically indicated – up to 7 days use
	10 x 10cm	5		
	12.5 x 12.5cm	5		
ActivHeal Foam adhesive	7 x 7cm	5		<b>Only if Aquacel Foam Pro adhesive® not suitable</b> If clinically indicated – up to 7 days use
	10 x 10cm	5		
	12.5 x 12.5cm	5		
Aquacel foam non adhesive	10 x 10cm	5		If clinically indicated, up to 7 days use
	15 x 15cm	5		
Duoderm® extra thin	10 x 10cm	5		<b>DO NOT apply to diabetic/ischaemic foot ulcers Avoid in infected wounds</b> Up to 7 days use
Aquacel® Extra	5 x 5cm	5		If clinically indicated – up to 7 days
	10 x 10cm	5		
	15 x 15cm	5		
Aquacel® ribbon	2 x 45cm	2		If clinically indicated – up to 7 days
Cosmopore® E	5 x 7.2 cm	10		If clinically indicated – up to 7 days
	8 x 10 cm	10		
Atrauman®	7.5 x 10cm	5		Effective for up to 7 days
	10 x 20cm	5		
Derma-Protective® Plus Skin Protectant	115g	1		No specific guidance regarding duration – use clinical judgement
365 film® dressings	6 x 7cm	10		Dressings can be left in place for up to 7 days, if signs of infection are observed, remove the dressing.
	10 x 12cm	10		
Tubifast®	Blue 7.5cm x 3m	5		No specific guidance regarding duration– use clinical judgement
	Blue 7.5cm x 5m	5		
	Yellow 10.75cm x 3m	5		
	Yellow 10.75cm x5m	5		
K- soft®	10cm x 3.5 m	5		No specific guidance regarding duration– use clinical judgement

K-Soft® Long	10cm x 4.5m	5		No specific guidance regarding duration– use clinical judgement
K- lite®	10cm x 4.5 m	5		No specific guidance regarding duration – use clinical judgement
K-Plus®	10 x 8.7cm	5		No specific guidance regarding duration – use clinical judgement
Kerrapro® sheet	10 x 10 x 1.2	1		Reusable - can be washed with soap and water, dried thoroughly and reused on the same patient.

Please indicate the dressings that are no longer required:

\*Additional information and maximum number of days dressings can stay on are for general reference only. Please exercise holistic clinical judgement when prescribing

***\*\*For full list of formulary items, please see the Greenwich Combined Wound Care Formulary***